

# COMPARE-1 CLINICAL TRIAL

PROTOCOL PRE-SPECIFIED INTERIM ANALYSIS OF FIRST 150 PATIENTS  
(COMPARE PILOT) AFTER 12-MONTHS OF FOLLOW-UP PRESENTED AT LINC 2018<sup>1</sup>

## TRIAL DESIGN

**WORLD'S FIRST** Head-to-Head Prospective, **RANDOMIZED (1:1) CONTROLLED**

Trial comparing the lower dose Paclitaxel RANGER DCB (2 µg/mm<sup>2</sup>) to the higher dose Paclitaxel IN.PACT DCB (3.5 µg/mm<sup>2</sup>)

## HYPOTHESIS

1. Lower dose Paclitaxel Ranger DCB (2 µg/mm<sup>2</sup>) with TransPax excipient is able to achieve similar patency as compared to the higher dose Paclitaxel IN.PACT DCB (3.5 µg/mm<sup>2</sup>) with Urea excipient
2. Post-procedure outcomes with the 0.018" Ranger DCB platform are similar to outcomes with the 0.035" and 0.018" IN.PACT DCB platform

- **No statistically significant difference in patency between the higher dose Paclitaxel IN.PACT DCB and the lower dose Paclitaxel Ranger DCB (89% and 84%, respectively, p = 0.7).**

- **Ranger is loaded with 2 µg/mm<sup>2</sup> and IN.PACT is loaded with 3.5 µg/mm<sup>2</sup> of Paclitaxel**

BASELINE CHARACTERISTICS AND RESULTS	Ranger™ (n=74)	IN.PACT™ (n=76)	p-value
Total Occlusions	<b>39%</b>	<b>45%</b>	0.5
Total Occlusion Length	<b>111 mm</b>	<b>95 mm</b>	0.5
Target Lesion Length	<b>117 mm</b>	<b>122 mm</b>	0.8
Moderate to Severe Calcification	<b>58%</b>	<b>61%</b>	0.7
Diabetics	<b>34%</b>	<b>37%</b>	0.7

<sup>1</sup>Results from the 150 patients from the pilot phase. Overall trial will enroll up to 414 patients.

\*KM Estimate

# COMPARE-1 CLINICAL TRIAL

## TRIAL DETAILS

### OBJECTIVE:

To compare two different Paclitaxel coated balloons (with different coatings and different Paclitaxel dose density) in the treatment of high grade stenotic or occluded lesions in SFA and/or PPA in PAD patients with Rutherford class 2-4.

### TRIAL DESIGN:

- Prospective, multicenter, randomized trial
- Randomization 1:1
- Phase 1: Pilot Study (150 patients)
- Phase 2: Extension (up to 414 patients) for testing of a formal non-inferiority hypothesis
- Stratification according to lesion length
- Follow-up clinical visits at 6, 12, 24 months

	Ranger™	IN.PACT™
Drug Dose	2.0 µg/mm <sup>2</sup>	3.5 µg/mm <sup>2</sup>
Longest Length Available	200 mm	150 mm
Platform	Sterling™	Admiral™ or Pacific™
Excipient	TransPax™	Urea

### KEY ENROLLMENT CRITERIA:

- Rutherford 2, 3 or 4
- Stenotic, restenotic or occlusive lesions (≥70% stenosis) in the native non-stented SFA/PPA
- No prior treatment with drug coated balloons or drug-eluting stents in the treated limb
- Lesion ≤300 mm, RVD ≥ 4 mm and ≤ 6.5 mm

### PRIMARY ENDPOINT DEFINITIONS:

#### Patency — Core lab adjudicated (12 months):

Efficacy: patency rate after one year defined as absence of clinically driven Target Lesion Revascularization (due to symptoms and drop of ABI of ≥ 20% or > 0.15 when compared to post-procedure) or restenosis with Peak Velocity Ratio > 2.4 evaluated by Duplex Ultrasound

**Boston  
Scientific**  
Advancing science for life™

**Peripheral Interventions**  
300 Boston Scientific Way  
Marlborough, MA 01752-1234  
[www.bostonscientific.com](http://www.bostonscientific.com)

To order product or for more information  
contact customer service at 1.888.272.1001.

© 2018 Boston Scientific Corporation  
or its affiliates. All rights reserved.

PI-529803-AA MAR2018

**CAUTION:** Ranger is an Investigational Device. Limited by Federal (or US) law to investigational use only.  
Not available for sale.

**Ranger™** Paclitaxel Coated PTA Balloon Catheter is manufactured by Hemoteq AG. All cited trademarks are the property of their respective owners.